



CARE SURGICAL

# CS Comfort Shin Supports (Pair)

## Instructions for use

- + Ensure feet do not touch operating table
- + Ankles should extend beyond Shin support



Manufacturer

UK: Care Surgical Ltd  
 Unit 6 Ringtail Road  
 Burscough, L40 8JY  
 Tel: 01704 336671

USA: Care Surgical LLC  
 Marlborough  
 MA 01752  
 Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)

Care Surgical Supplies Limited, 69 Esker Woods Drive  
 Lucan, Co Dublin, K78PX45, Ireland

EC REP



Medical Device



Reuse



CE Mark



Latex Free



Non Sterile



5 060624 120725

Australian Sponsor, Summit Surgical  
 PO Box 60, Landsborough  
 QLD 4550 Australia

LOT Number:

Product Code: CSM-2820 Rev 06 08/21

# Head-to-Toe Check list

This check should be completed once the patient is in the prone position and before the surgical drapes are applied. At all stages, check that no tubing is trapped between the patient and the supports, and that all bony protrusions are padded.

## The head

- + Eyes lightly taped, no padding or goggles
- + Face supported on the cheeks and forehead
- + No pressure on the nerves around the eyes
- + No pressure or minimal pressure on the chin




## Neck

- + Chin and throat not touching the chest support
- + Natural neutral neckline

## Chest, abdomen and pelvis

- + ECG electrodes positioned on the patient's back
- + Breasts positioned medially with minimum pressure
- + Adequate space to avoid compression of the abdomen
- + Male genitalia free of the pelvic support



## The arms

- + Arms forward and not abducted beyond 90°, or arms secure by the patient's side
- + No pressure on the axillary nerve
- + No body weight supported on the arms

To see the full range of Care Surgical Medical Products, visit [www.care-surgical.com](http://www.care-surgical.com)



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
Tel: 01704 336671

USA: Care Surgical LLC  
Marlborough  
MA 01752  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)

Care Surgical Supplies Limited, 69 Esker Woods Drive  
Lucan, Co Dublin, K78PX45, Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia

Product Code: CSM-2820 Rev 06 08/21



## GTIN

- + 5060624120725

## Description

- + Shin supports made from polyurethane foam with a soft wipeable top and non-slip base

## Intended use

- + Support patients' shins in the prone position

## Indications for Use

- + Support the patient's knees in the prone position
- + Reduce the risk of pressure ulcers
- + Latex and phthalate free
- + Reusable: Clean using non-alcohol hospital wipes
- + Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)
- + Instructions for use are packaged with the product

## Warnings, Precautions, Contraindications & Adverse Effects

- + Incorrect positioning of the patient on this product could result in the patient experiencing:
  - + Pressure ulcers
  - + Neurological damage, which can lead to temporary or permanent paralysis
  - + Venous pressure, which can lead to ischaemia
- + Examine the product for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced
- + Do not immerse in water
- + Do not autoclave

## Packaging

- + Each Unit is packaged in a recyclable polyethene bag
- + 1 pair is packaged into a cardboard box for sale

## Transport & Storage

- + Transport and store in dry conditions and out of direct light where possible
- + Avoid storing where sharp objects could pierce the skin of the product