



Technical Data Sheet

CSM2908 CS Plexus Wedges (Pair)



Description

- Reusable
- Placed under the patients shoulder blades in supine position the wedges help to reduce the pressure on the brachial plexus and to improve the comfort
- Comfortable for both conscious and anaesthetised patients

Application Area

- OR
- Anaesthesia

Measurements

- Length: 195mm/7.68" (Per Wedge)
- Width: 95mm/3.74" (Per Wedge)
- Height: 10mm - 70mm (Per Wedge, lowest height to tallest height)

Material Composition

- Polyurethane foam core
- PU Cover
- Latex & phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2908	CS Plexus Wedges (Pair)	1 Pair

Packaging Information

- Each Pair is packaged in a recyclable polyethene bag
- 1 pair is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can lead to temporary or permanent paralysis
- Venous pressure, which can lead to ischaemia. Examine the product for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624125966

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark