Technical Data Sheet

CSM3043 CS Fixed Blade Clamp



Order Information

Product Code	Product Name	Qty
CSM3043	CS Fixed Blade Clamp	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Description

- For use on compatible surgical tables in surgical environments and ICUs
- Reusable - clean using alcohol-based hospital wipes
- Latex and phthalate free
- Intended to fix or mount blade-style instrumentation to operating table rails to the device
- Suitable for 25*10mm (EU), 28.5*9.5mm (US), and 30*8mm (JPN) side rails

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Transport and store in dry conditions

Application Area

- OR

Contraindications

- Incorrect attachment to the surgical tables could result in the attachments failing
- Do not use if the blade clamp appears to be defective or damaged
- Only to be be used by trained medical personnel
- Do not exceed the device safe working load
- Do not immerse in water
- Do not autoclave
- After each use, clean the device with alcohol-based wipes
- Make sure that the device is dry before you store it or use it again

Measurements

- Length 75mm/2.95"
- Width 45mm/1.77"
- Height 50mm/1.97"

Material Composition

- Main Body: Stainless Steel
- Screw: Stainless Steel & Plastic (ABS)
- Latex & phthalate free

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624126697

Labelling

Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com

EU REP

Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Medical Device | Reuse | Latex Free | Non Sterile | UK Conformity Assessed | CE Mark