




CARE SURGICAL

CS Disposable Lateral Head Support (Adult)



It is recommended that the patient's head is turned every 2 to 4 hours per ICS guidelines

**No need to disconnect the tubing
Single use to aid infection control**



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
Marlborough
MA 01752
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, 69 Esker Woods Drive
Lucan, Co Dublin, K78PX45, Ireland



Do not reuse



CE Mark



Latex Free



Non Sterile



Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

LOT Number:

Product Code: CSM-2410 Rev 4 06/21

GTIN

- + 5060624120800

Description

- + Horseshoe-shaped head support made from polyurethane foam

Intended use

- + Supports patient's heads in the lateral position

Indications for Use

- + Single use - disposal after use with regular clinical waste
- + Recommended for adults, but use will depend on patient's stature
- + Supports patient's head in the lateral position
- + Can be used for patients who are in the prone or lateral position
- + Reduces the incidence of pressure ulcers and marking
- + Patient can be repositioned without disconnecting their tubing, so PEEP pressure is maintained
- + Latex and phthalate free
- + Instructions for use are packaged with the product

Warnings, Precautions, Contraindications & Adverse Effects

- + There is a 1/11,000 risk of permanent visual loss associated with prone positioning, irrespective of the head support used
- + Incorrect positioning of the patient in the this product can increase the risk of:
 - + POVL
 - + Pressure ulcers
 - + Neurological damage
 - + Venous pressure, which can lead to ischaemia
- + Do not use if this product is torn or damaged

Packaging

- + Each Unit is packaged in a recyclable polyethene bag
- + 20 units are packaged into a cardboard box for sale

Transport & Storage

- + Transport and store in dry conditions and out of direct light where possible
- + Exposure to light may cause yellowing of areas of the product. This is purely cosmetic and will not affect the quality or the safety of the product