



# Technical Data Sheet

## CSM-2510 CS Prone Face Cushion



### Description

- Prone Face Cushion - Made from polyurethane foam
- For use with the CS Head Support System and other compatible head support systems
- Distributes pressure on the patient's cheeks and forehead while in the prone position
- For patients aged 8 and over, dependent on stature
- For patients weighing up to 200kg

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 312mm/12.28"
- Width 224mm/8.82"
- Height 146mm/5.75"

### Material Composition

- Polyurethane foam
- Latex & phthalate free

### Order Information

| Product Code | Product Name          | Qty      |
|--------------|-----------------------|----------|
| CSM-2510     | CS Prone Face Cushion | 10 Units |

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 10 units are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- POVL
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120015

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson,  
MA 01749, USA  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main  
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity  
Assessed



CE Mark

Rev.02 - 08/2024