




# CARE SURGICAL

# CS Comfort Covers



## CS Comfort Covers - 4 Hip/Thigh Covers & 1 Chest Cover

Enhances infection control and protects the pads from damage. Reduces the chance of pressure related injury.



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
Tel: 01704 336671

USA: Care Surgical LLC  
Marlborough  
MA 01752  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)

Care Surgical Supplies Limited, 69 Esker Woods Drive  
Lucan, Co Dublin, K78PX45, Ireland

EC REP



Do not reuse



CE Mark



Latex Free



Non Sterile



5 060624 120220

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia

LOT Number:

Product Code: CSM-2600-A Rev 07 08/21

## GTIN

- + 5060624120916

## Description

- + Disposable Drapes with waterproof lining and 30mm foam top

## Intended use

- + Protective covers

## Indications for Use

- + Single Use
- + For use on the CS Prone Patient Supports or Allen Advance Table supports
- + Helps to protect the chest, hip and thigh pads on the Allen Advance table from wear and tear
- + Aids infection control
- + 30mm foam padding to reduce the risk of the patient experiencing marking or pressure ulcers
- + Each set contains 1x chest, 2x hip and 2x thigh covers
- + Latex and phthalate free
- + Instructions for use are packaged with the product

## Warnings, Precautions, Contraindications & Adverse Effects

- + If incorrectly applied to the pads, this product can cause pressure marks or ulcers on the patient

## Packaging

- + Each Unit is packaged in a recyclable polyethene bag
- + 6 sets are packaged into a cardboard box for sale

## Transport & Storage

- + Transport and store in dry conditions and out of direct light where possible
- + Exposure to light may cause yellowing of areas of the product. This is purely cosmetic and will not affect the quality or the safety of the product