




CARE SURGICAL

CS Jackson Comfort Covers 30mm 5-piece



CS Comfort Covers - 4 hip/Thigh Covers & 1 Chest Cover Compatible with Jackson Table

Enhances infection control and protects the pads from damage. Reduces the chance of pressure related injury.



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
Marlborough
MA 01752
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, 69 Esker Woods Drive
Lucan, Co Dublin, K78PX45, Ireland

EC REP



Do not reuse



CE Mark



Latex Free



Non Sterile



5 060624 120350

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

LOT Number:

Product Code: CSM-2602 Rev 02 06/21

GTIN

- + 5060624120350

Description

- + Disposable Drapes with waterproof lining

Intended use

- + Protective covers

Indications for Use

- + Single Use
- + For use on the Jackson Table
- + Helps to protect the Jackson Table pads from wear and tear
- + Aids infection control
- + Softly padded chest pad and lightly padded hip and thigh covers reduce the risk of the patient experiencing marking or pressure ulcers

Set consists of:

- + A soft padded cover for the chest pad
- + 4 hip & Thigh covers for the hip and thigh pads
- + Latex and phthalate free
- + Instructions for use are packaged with the product

Warnings, Precautions, Contraindications & Adverse Effects

- + If incorrectly applied to the pads of the Jackson Table, this product can cause pressure marks or ulcers on the patient

Packaging

- + Each Unit is packaged in a recyclable polyethene bag
- + 6 sets are packaged into a cardboard box for sale

Transport & Storage

- + Transport and store in dry conditions and out of direct light where possible
- + Exposure to light may cause yellowing of areas of the product. This is purely cosmetic and will not affect the quality or the safety of the product