




CARE SURGICAL

CS Comfort Covers Large 30mm - 5 Piece



CS Comfort Covers Large
2 Large Hip Covers
2 Thigh Covers
1 Chest Cover

Enhances infection control and protects the pads from damage. Reduces the chance of pressure related injury.



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
Marlborough
MA 01752
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, 69 Esker Woods Drive
Lucan, Co Dublin, K78PX45, Ireland

EC REP



Do not reuse



CE Mark



Latex Free



Non Sterile

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

LOT Number:

Product Code: CSM-2603 Rev 02 06/21

GTIN

- + 5060624121203

Description

- + Disposable Drapes with waterproof lining and 30mm foam top

Intended use

- + Protective covers

Indications for Use

- + Single use - disposal after use with regular clinical waste
- + For use on a number of modular tables and frames including the CS Spine Frame and the Allen® Advance table
- + Includes 1 chest cover 2 hip covers and 2 thigh covers
- + Waterproof backing helps to protect the padded supports of the table/frame from wear and tear
- + Aids infection control
- + 30mm foam padding to reduce the risk of the patient experiencing marking or pressure ulcers.
- + Latex and phthalate free
- + Instructions for use are packaged with the product

Warnings, Precautions, Contraindications & Adverse Effects

- + Incorrect positioning of the cover on the chest support, or of the patient on this product can cause:
 - + Pressure marking or ulcers
 - + Neurological damage, which can cause paraplegia;
 - + Venous pressure, which can lead to ischaemia

Packaging

- + Each Unit is packaged in a recyclable polyethene bag
- + 6 sets are packaged into a cardboard box for sale

Transport & Storage

- + Transport and store in dry conditions and out of direct light where possible
- + Exposure to light may cause yellowing of areas of the product. This is purely cosmetic and will not affect the quality or the safety of the product