




CARE SURGICAL

CS Jackson Standard Covers - 3 piece



**2 Hip/Thigh Covers & 1 Chest
Cover - Compatible with
Jackson style supports**

Single use to aid infection control.
Protects the pads from damage.
Reduces the chance of pressure
related injury.



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
Marlborough
MA 01752
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, 69 Esker Woods Drive
Lucan, Co Dublin, K78PX45, Ireland

EC REP



Do not reuse



CE Mark



Latex Free



Non Sterile



Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

LOT Number:

Product Code: **CSM-2604** Rev 02 06/21

GTIN

- + 5060624120398

Description

- + Disposable Drapes

Intended use

- + Protective covers

Indications for Use

- + Single Use
- + For use on the Jackson Table
- + Helps to protect the Jackson Table pads from wear and tear
- + Aids infection control
- + Softly padded chest pad and lightly padded hip and thigh covers reduce the risk of the patient experiencing marking or pressure ulcers

Set consists of:

- + A padded cover for the chest pad
- + 2 hip & Thigh covers for the hip and thigh pads
- + Latex and phthalate free
- + Instructions for use are packaged with the product

Warnings, Precautions, Contraindications & Adverse Effects

- + If incorrectly applied to the pads of the Jackson Table, this product can cause pressure marks or ulcers on the patient

Packaging

- + Each Unit is packaged in a recyclable polyethene bag
- + 6 sets are packaged into a cardboard box for sale

Transport & Storage

- + Transport and store in dry conditions and out of direct light where possible
- + Exposure to light may cause yellowing of areas of the product. This is purely cosmetic and will not affect the quality or the safety of the product