




CARE SURGICAL

CS Comfort Covers - 30mm 3-Piece



CS Comfort Covers - 2 Hip/Thigh Covers & 1 Chest Cover

Enhances infection control and protects the pads from damage. Reduces the chance of pressure related injury.



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
Marlborough
MA 01752
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, 69 Esker Woods Drive
Lucan, Co Dublin, K78PX45, Ireland



Do not reuse



CE Mark



Latex Free



Non Sterile



Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

LOT Number:

Product Code: CSM-2606-1 Rev 09 06/21

GTIN

- + 5060624121166

Description

- + Disposable Drapes with waterproof lining and 30mm foam top

Intended use

- + Protective covers

Indications for Use

- + Single Use
- + For use on the CS Prone Patient Supports or Allen Advance Table supports
- + Helps to protect the chest and hip pads on the Allen Advance table from wear and tear
- + Aids infection control
- + 30mm foam padding to reduce the risk of the patient experiencing marking or pressure ulcers
- + Each set contains 1x chest and 2x hip covers
- + Latex and phthalate free
- + Instructions for use are packaged with the product

Warnings, Precautions, Contraindications & Adverse Effects

- + If incorrectly applied to the pads, this product can cause pressure marks or ulcers on the patient

Packaging

- + Each Unit is packaged in a recyclable polyethene bag
- + 6 sets are packaged into a cardboard box for sale

Transport & Storage

- + Transport and store in dry conditions and out of direct light where possible
- + Exposure to light may cause yellowing of areas of the product. This is purely cosmetic and will not affect the quality or the safety of the product