



Technical Data Sheet

CSM-2606-1 CS Comfort Covers - 3 Piece



Description

- Disposable Foam Drapes with waterproof lining
- For use on the CS Spine Frame or Allen Advance table
- Foam pads and water-resistant drapes help to protect the chest and hip pads on the patient supports from wear and tear
- Single-use aids infection control
- Light padding to reduce the risk of the patient experiencing marking or pressure ulcers
- Foam covering absorbs excess fluids
- Water-resistant drape protects the pads
- Loose-fitting drape moves with the patient reducing friction and helping to prevent skinshear

Application Area

- OR
- Anaesthesia

Measurements

Chest

- Length 480mm/18.89"
- Width 275mm/10.82"
- Height 10mm/0.39"

Hips

- Length 207mm/8.14"
- Width 162mm/6.37"
- Height 10mm/0.39"

*10mm foam pad with water resistant drape attached

Order Information

Product Code	Product Name	Qty
CSM-2606-1	CS Comfort Covers - 3 Piece	6 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- If incorrectly applied to the pads, this product can cause pressure marks or ulcers on the patient
- Do not reuse
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122309

Material Composition

- White foam: Polyurethane
- Blue non-woven (waterproof): 55% Woodpulp 45% Polyester
- Latex & phthalate free

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 16/10/2024