



**CARE SURGICAL**  
Specialists in Patient Positioning

# CS Head Tray Attachment

The CS Head Tray Attachment - Lowered is compatible with most standard dual-rail operating systems on the market.

- + Radiolucent design
- + Quick-release clipping mechanism for easy attachment and adjustment



\*Compatible with rail sizes that measure 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall with 14.5" (36.8cm) space between.

## CS Head Tray Attachment

Code: CSM-2714  
Qty: 1 Unit

## Related Products:



**CS Prone Plus Head Support System**

Code: CSM-2565  
Qty: 1 System  
(Helmet & Mirror)



**CS Prone Head Support System**

Code: CSM-2560  
Qty: 1 System  
(Helmet & Mirror)




**CS Prone Plus Face Cushion**

Code: CSM-2512  
Qty: 10 Units



**CS Prone Small Face Cushion**

Code: CSM-2514  
Qty: 10 Units



**CS Prone Face Block 7" - Left & Right Intubation Slot (compression packed)**

Code: CSM-2536V  
Qty: 10 Units



**Contact us for more information:**  
[www.care-surgical.com](http://www.care-surgical.com) | [info@care-surgical.com](mailto:info@care-surgical.com)