



**CARE SURGICAL**  
Specialists in Patient Positioning

# CS Prone Chest and Pelvic Rolls

Raises the patient's chest and pelvis in the prone position to allow free movement of the abdomen.



## CS Prone Chest and Pelvic Rolls

Code: CSM-2902

Qty: 2

- + Assist to achieve neutral neckline
- + Soft polyurethane foam
- + Single-use

## Related Products:




**CS ProneOne Head Support**

Code: CSM-2520  
Qty: 8 Units



**CS ProneOne Head Support - Small**

Code: CSM-2519  
Qty: 10 Units



**CS Prone Face Block 7\"/>**

Code: CSM-2536V  
Qty: 10 Units



**CS Standard Arm Cushions**

Code: CSM-2575  
Qty: 10 Units



**CS Adjustable Prone Leg Cushion**

Code: CSM-2590  
Qty: 2 Unit



**Contact us for more information:**  
[www.care-surgical.com](http://www.care-surgical.com) | [info@care-surgical.com](mailto:info@care-surgical.com)