




CARE SURGICAL

CS ProAxis® Compatible Kit



1 Chest Cushion &
2 Hip Hinge Covers
Compatible with
ProAxis® Table

Enhances infection control and protects
the pads from damage. Reduces the
chance of pressure related injury.



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
Marlborough
MA 01752
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, 69 Esker Woods Drive
Lucan, Co Dublin, K78PX45, Ireland

EC REP



Do not reuse



CE Mark



Latex Free



Non Sterile



5 060624 120473

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

LOT Number:

Product Code: CSM-4570 Rev 03 06/21

GTIN

- + 5060624120473

Description

- + Disposable, waterproof drapes

Intended use

- + Protective covers

Indications for Use

- + Single Use
- + For use on the OSI ProAxis Table
- + Helps to protect the chest and hip pads on the ProAxis Table from wear and tear
- + Aids infection control
- + Light padding on the chest piece to reduce the risk of the patient experiencing marking or pressure ulcers.
- + Latex and phthalate free

Warnings, Precautions, Contraindications & Adverse Effects

- + If incorrectly applied to the pads, the CS ProAxis Table Compatible Covers can cause pressure marks or ulcers on the patient

Packaging

- + Each Unit is packaged in a recyclable polyethylene bag
- + 6 sets are packaged into a cardboard box for sale

Transport & Storage

- + Transport and store in dry conditions and out of direct light where possible
- + Exposure to light may cause yellowing of areas of the product. This is purely cosmetic and will not affect the quality or the safety of the product