



Technical Data Sheet

CSM2570 CS Bow & Wilson Standard Covers



Description

- Disposable Drapes
- For use on the Wilson or Bow Frame
- Helps to protect the Wilson or Bow Frame from wear and tear
- Single us to aid infection control
- Kit consists of:
 - Sleeves to cover the frame's bar
 - Square drape to cover the frame's base
 - Bar roll to cover the frame's central bar

Application Area

- OR
- Anaesthesia

Measurements

- Length 760mm/29.92"
- Width 250mm/9.84"
- Height 5mm/0.20"

Material Composition

- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2570	CS Bow & Wilson Standard Covers	5 Units

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 5 sets are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- POVL
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Only to be used by trained medical personnel
- Do not use if this product is torn or damaged

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120176

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark