



CARE SURGICAL

Specialists in Patient Positioning

INSTRUCTIONS FOR USE

CSM2821 **CS Knee Rings (Paar)**




CARE SURGICAL

Product Code: CSM2821

CS Knee Rings (Pair)

Instructions for use

- + Place the knees in the centre of the rings for minimum pressure



⚠ Only to be used by trained clinicians.

Ensure product is secure before positioning patient.


Cleaning instructions:



Clean using standard hospital wipes



Do not immerse in water



Do not autoclave



UK: Care Surgical Ltd
Unit 6 Ringtail Road,
Burscough,
L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



SAFE DISPOSAL: customers should adhere to any local laws and regulations to the safe disposal of medical devices.
DO NOT USE products that appear to be damaged or worn.
In an event of product damage or patient injury, please contact Care Surgical.



5 060624 121128

EU REP




CARE SURGICAL

Artikelnummer: CSM2821

CS Knee Rings (Paar)

Gebrauchsanweisung

- + Platzieren Sie die Knie in der Mitte der Ringe für minimalen Druck



⚠ Nur für den Einsatz durch geschultes Personal geeignet
Vergewissern Sie sich, dass das Produkt richtig liegt,
bevor Sie den Patienten positionieren.


Reinigungshinweise:



Mit handelsüblichen Kranken-
haus-Reinigungstüchern
reinigen



Nicht in Wasser eintauchen



Nicht autoklavieren



UK: Care Surgical Ltd
Unit 6 Ringtail Road,
Burscough,
L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



Maximales
Patientengewicht



SICHERE ENTSORGUNG: Kunden sollten alle lokalen Gesetze und
Vorschriften zur sicheren Entsorgung von Medizinprodukten einhalten.
VERWENDEN SIE KEINE PRODUKTE die beschädigt oder abgenutzt
erscheinen. Im Falle einer Produktbeschädigung oder einer Verletzung
des Patienten wenden Sie sich bitte an Care Surgical.



5 060624 121128

EU REP



CARE SURGICAL
Specialists in Patient Positioning