




Technical Data Sheet

CSM2536 CS Prone Face Block 7" -Left & Right Intubation Slot



Description

- Supports patients' faces in the prone position
- Stand alone/ all-in-one product
- Firm base provides stability
- Soft foam top for patient comfort and protection
- Distributes the pressure on the patient's cheeks, forehead and chin while in the prone position
- Reduces the incidence of pressure ulcers and marking
- Raises the patient's head to assist with achieving a natural, neutral neckline
- Hole in the side of the block to hold tubing
- For patients aged 8 and over, dependent on stature
- For patients weighing up to 200kg

Application Area

- OR
- Anaesthesia

Measurements

- L270 | W246 | H178mm
- L10.6 | W9.6 | H7"

Material Composition

- Made from a layer of soft polyurethane foam & firm polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2536	CS Prone Face Block 7" -Left & Right Intubation Slot	8 Units

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 8 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- POVL
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Only to be used by trained medical personnel
- Do not use if this product is torn or damaged

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120077

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

Rev.02 01/04/2025