



Technical Data Sheet

CSM-2410 CS Disposable Lateral Head Support (Adult)



Description

- Horseshoe-shaped head support made from polyurethane foam
- Recommended for adults, but use will depend on patient's stature
- Supports patient's head in the lateral position
- Can be used for patients who are in the supine or lateral position
- Reduces the incidence of pressure ulcers and marking
- Patient can be repositioned without disconnecting their tubing, so PEEP pressure is maintained

Application Area

- OR
- Anaesthesia

Measurements

- Length 200mm/7.87"
- Width 200mm/7.87"
- Height 70mm/2.76"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2410	CS Disposable Lateral Head Support (Adult)	20 units

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 20 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient in the this product can increase the risk of:

- POVL
 - pressure ulcers
 - neurological damage
 - venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
 - Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120800

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023