



Technical Data Sheet

CSM-2512 CS Prone Plus Face Cushion



Description

- Cushion to position the patients face in prone
- universal face cut-out (suitable for patients from 8 years upwards)
- Optimal pressure distribution on forehead and cheeks
- No pressure on the sensitive chin area
- Raises the patient's head to assist with achieving a natural, neutral neckline
- For use in combination with the CS Prone Plus Head Support System (CSM-2565, consisting of helmet and mirror) or other compatible head support systems
- Free of latex

Application Area

- OR
- Anaesthesia
- ICU

Measurements

- L312 | W224 | H146 mm
- L12.28 | W8.82 | H5.75"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2512	CS Prone Plus Face Cushion	10 units

Packaging Information

- Each Unit is packaged in a recyclable polyethylene bag
- 10 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120008

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.02 10/23