



Technical Data Sheet

CSM2530 CS Prone Face Mask



Description

- Supports patients' faces in the prone position
- Compatible with the Allen C-Prone Head Positioner
- Distributes the pressure on the patient's cheeks, forehead and chin whilst in the prone position
- Reduces the incidence of pressure ulcers and marking

Application Area

- OR
- Anaesthesia

Measurements

- Length: 330mm/12.99"
- Width: 280mm/11.02"
- Height: 50mm/1.97"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2530	CS Prone Face Mask	5 Units

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 5 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- POVL
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624120572

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark