




# Technical Data Sheet

## CSM-2536V CS Prone Face Block 7" -Left & Right Intubation Slot (Compression Packed)



### Description

- Supports patients' faces in the prone position
- Stand alone/all-in-one product
- Firm base provides stability
- Soft foam top for patient comfort and protection
- Distributes the pressure on the patient's cheeks, forehead and chin while in the prone position
- Reduces the incidence of pressure ulcers and marking
- Raises the patient's head to assist with achieving a natural, neutral neckline
- Hole in the side of the block to hold tubing
- For patients aged 8 and over, dependent on stature
- For patients weighing up to 200kg

### Application Area

- OR
- Anaesthesia

### Measurements

- L270 | W246 | H178mm
- L10.6 | W9.6 | H7"

### Material Composition

- Made from a layer of soft polyurethane foam & firm polyurethane foam
- Latex & phthalate free

### Order Information

Product Code	Product Name	Qty
CSM-2536V	CS Prone Face Block 7" -Left & Right Intubation Slot (Compression Packed)	10 Units

### Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 10 units are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- POVL
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Only to be used by trained medical personnel
- Do not use if this product is torn or damaged

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120732

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023