




Technical Data Sheet

CSM-2551V CS Prone Face Block 5" - Left & Right Intubation Slot (Compression Packed)



Description

- Supports patients' faces in the prone position
- Stand alone/ all-in-one product
- Firm base provides stability
- Memory foam top for patient comfort and protection
- Distributes the pressure on the patient's cheeks, forehead and chin while in the prone position
- Reduces the incidence of pressure ulcers and marking
- Raises the patient's head to assist with achieving a natural neutral neckline
- Hole in the side of the block to hold tubing
- For patients aged 8 and over, dependent on stature
- For patients weighing up to 200kg
- Remove from bag 1 (one) hour before use to allow the cushion to regain required shape and size

Application Area

- OR
- Anaesthesia

Measurements

- Length 270mm/10.63"
- Width 246mm/9.69"
- Height 127mm/5.00"

Material Composition

- Made from a layer of memory foam & polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2550	CSM-2551V CS Prone Face Block 5" - Left & Right Intubation Slot (Compression Packed)	14 Units

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 14 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia;
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120091

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023