



Technical Data Sheet

CSM-2570 CS Bow & Wilson Standard Covers



Description

- Disposable Drapes
- For use on the Wilson or Bow Frame
- Helps to protect the Wilson or Bow Frame from wear and tear
- Single us to aid infection control
- Kit consists of:
 - Sleeves to cover the frame's bar
 - Square drape to cover the frame's base
 - Bar roll to cover the frame's central bar

Application Area

- OR
- Anaesthesia

Measurements

- Length 760mm/29.92"
- Width 250mm/9.84"
- Height 5mm/0.20"

Material Composition

- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2570	CS Bow & Wilson Standard Covers	5 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 5 sets are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight

Contraindications

- If incorrectly applied to the Frame, this product can cause pressure marks or ulcers on the patient
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120176

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not
Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark

Rev.01 11/09/2023