




# Technical Data Sheet

## CSM-2585 CS Egg Box Foam 80" x 20'



### Description

- Padding for patients
- Can be used on operating tables, trolleys beds or chairs
- Each sheet is 200x50cm
- Can be cut to size, or tear along perforations
- Reduces and distributes pressure to reduce the risk of pressure ulcers.
- The undulating structure of the foam enables air to circulate to assist in maintaining the patient's temperature.

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 2032mm/80"
- Width 508mm/20"
- Height 15mm/0.59"

### Material Composition

- Polyurethane foam core
- Latex & phthalate free

### Order Information

Product Code	Product Name	Qty
CSM-2585	CS Egg Box Foam 80" x 20'	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

- Incorrect placement of this product could cause pressure marks or ulcers on the patient
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120190

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023