



Technical Data Sheet

CSM2600 CS Comfort Covers



Description

- Disposable Foam Drapes with waterproof lining
- For use on the CS Prone Patient Supports or Allen Advance
- Foam pads and water-resistant drapes help to protect the chest, hip and thigh pads on the patient supports from wear and tear
- Single-use aids infection control
- Light padding to reduce the risk of the patient experiencing marking or pressure ulcers
- Foam covering absorbs excess fluids
- Water-resistant drape protects the pads
- Loose-fitting drape moves with the patient reducing friction and helping to prevent skinshear

Application Area

- OR
- Anaesthesia

Measurements

Chest

- Length: 480mm/18.89"
- Width: 275mm/10.82"
- Height: 10mm/0.39"*

Hips

- Length: 207mm/8.14"
- Width: 162mm/6.37"
- Height: 10mm/0.39"*

Thighs

- Length: 207mm/8.14"
- Width: 162mm/6.37"
- Height: 10mm/0.39"*

*10mm foam pad with water resistant drape attached

Order Information

| Product Code | Product Name | Qty |
|--------------|-------------------|--------|
| CSM2600 | CS Comfort Covers | 6 Sets |

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- If incorrectly applied to the pads, this product can cause pressure marks or ulcers on the patient
- Do not reuse
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624121166

Material Composition

- White foam: Polyurethane
- Blue non-woven (waterproof): 55% Woodpulp
45% Polyester
- Latex & phthalate free

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA s
Tel: +1-866-243-4107

www.care-surgical.com

EU REP

Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Ref: CSM2600 - TDS (Rev.04) 03/2026