



Technical Data Sheet

CSM-2603 CS Comfort Covers Large- 5 piece



Description

- Disposable Foam Drapes with waterproof lining
- For use on the CS Prone Patient Supports or Allen Advance
- Foam pads and water-resistant drapes help to protect the chest, hip and thigh pads on the patient supports from wear and tear
- Single-use aids infection control
- Light padding to reduce the risk of the patient experiencing marking or pressure ulcers
- Foam covering absorbs excess fluids
- Water-resistant drape protects the pads
- Loose-fitting drape moves with the patient reducing friction and helping to prevent skinshear

Application Area

- OR
- Anaesthesia

Measurements

Chest:

- Length 480mm/18.89"
- Width 275mm/10.82"
- Height 10mm/0.39"

Hips:

- Length 310mm/12.20"
- Width 170mm/6.69"
- Height 10mm/0.39"

Thighs:

- Length 310mm/12.20"
- Width 170mm/6.69"
- Height 10mm/0.39"

* 10mm foam and water resistant drape attached

Order Information

Product Code	Product Name	Qty
CSM-2603	CS Comfort Covers Large- 5 piece	6 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624121203

Material Composition

- Waterproof lining
- 30mm foam top
- Latex & phthalate free

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023