




# Technical Data Sheet

## CSM2551V CS Prone Face Block 5" - Left & Right Intubation Slot (Compression Packed)



### Description

- Supports patients' faces in the prone position
- Stand alone/ all-in-one product
- Firm base provides stability
- Memory foam top for patient comfort and protection
- Distributes the pressure on the patient's cheeks, forehead and chin while in the prone position
- Reduces the incidence of pressure ulcers and marking
- Raises the patient's head to assist with achieving a natural neutral neckline
- Hole in the side of the block to hold tubing
- For patients aged 8 and over, dependent on stature
- For patients weighing up to 200kg
- Remove from bag 1 (one) hour before use to allow the cushion to regain required shape and size

### Application Area

- OR
- Anaesthesia

### Measurements

- Length: 270mm/10.63"
- Width: 246mm/9.69"
- Height: 127mm/5.00"

### Material Composition

- Made from a layer of memory foam & polyurethane foam
- Latex & phthalate free

### Order Information

Product Code	Product Name	Qty
CSM2551V	CSM-2551V CS Prone Face Block 5" - Left & Right Intubation Slot (Compression Packed)	14 Units

### Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 14 units are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia;
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624120091

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark