



# Technical Data Sheet

## CSM2565 CS Prone Plus Head Support System™



### Description

- Reusable
- Helmet holds the CS Prone Plus Face Cushion (CSM2512) or CS Small Face Cushion (CSM2514) to maintain their integrity throughout prone procedures
- Space for safe positioning of ET tubes.
- The helmet can be used with or without the mirror
- The mirror enables the clinician to monitor the patient's eyes and tubing during the procedure
- The cogs fitted to the mirror are used to adjust the height or angle of the helmet to achieve the safest natural neutral neckline

### Application Area

- OR
- Anaesthesia
- ICU

### Measurements

#### Helmet:

- Length: 317mm/12.48"
- Width: 233mm/9.17"
- Height: 155mm/ 6.10"

#### Mirror:

- Length: 220 mm/ 8.6"
- Width: 220 mm/8.6"

### Material Composition

- Helmet made from ABS Plastic
- Acrylic mirror, which incorporates 4 adjustable cogs
- Aluminium 6061

### Cleaning

- Clean with standard hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM2565	CS Prone Plus Head Support System™	1 System

### Packaging Information

- Each set is packaged in a recyclable polyethylene bag
- Each set is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- POVL
  - pressure ulcers
  - neurological damage
  - venous pressure, which can lead to ischaemia
- 
- Do not use if there is any damage on the helmet or mirror
  - Do not autoclave
  - Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Not fully radiolucent- cog screws may be visible on x-rays
- Sterility: N
- CE/FDA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624120121

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson,  
MA 01749, USA  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor:

MDSS AU Pty Ltd, Ground Floor, Suite 55,  
97-99 Bathurst Street, Sydney, NSW 2000, Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark