



Technical Data Sheet

CSM2602 CS Jackson Comfort Covers



Description

- Disposable Drapes with waterproof lining
- Helps to protect the Jackson Table pads from wear and tear
- Aids infection control
- Softly padded chest pad and lightly padded hip and thigh covers reduce the risk of the patient experiencing marking or pressure ulcers.
- Set consists of:
 - A soft padded cover for the chest pad
 - 4 hip & Thigh covers for the hip and thigh pads
 - Latex and phthalate free
 - Instructions for use are packaged with the product

Application Area

- OR
- Anaesthesia

Material Composition

- Waterproof lining
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2602	CS Jackson Comfort Covers	6 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- If incorrectly applied to the pads of the Jackson Table, this product can cause pressure marks or ulcers on the patient
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624120350

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark