



# Technical Data Sheet

## CSM2606-1 CS Comfort Covers - 3 Piece



### Description

- Disposable Foam Drapes with waterproof lining
- For use on the CS Spine Frame or Allen Advance table
- Foam pads and water-resistant drapes help to protect the chest and hip pads on the patient supports from wear and tear
- Single-use aids infection control
- Light padding to reduce the risk of the patient experiencing marking or pressure ulcers
- Foam covering absorbs excess fluids
- Water-resistant drape protects the pads
- Loose-fitting drape moves with the patient reducing friction and helping to prevent skinshear

### Application Area

- OR
- Anaesthesia

### Measurements

#### Chest

- Length: 480mm/18.89"
- Width: 275mm/10.82"
- Height: 10mm/0.39"

#### Hips

- Length: 207mm/8.14"
- Width: 162mm/6.37"
- Height: 10mm/0.39"

\*10mm foam pad with water resistant drape attached

### Order Information

Product Code	Product Name	Qty
CSM2606-1	CS Comfort Covers - 3 Piece	6 Sets

### Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

- If incorrectly applied to the pads, this product can cause pressure marks or ulcers on the patient
- Do not reuse
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624122309

### Material Composition

- White foam: Polyurethane
- Blue non-woven (waterproof): 55% Woodpulp 45% Polyester
- Latex & phthalate free

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark