



# Technical Data Sheet

## CSM2608 CS Bow & Wilson Comfort Covers



### Description

- For use on the Allen Bow or Wilson Frame
- Disposable Drape with 10mm foam and waterproof lining
- Helps to protect the Allen Bow or Wilson Frame from wear and tear
- Aids infection control
- Light padding to reduce the risk of the patient experiencing marking or pressure ulcers

### Application Area

- OR
- Anaesthesia

### Measurements

- Length: 760mm/29.92"
- Width: 250mm/9.84"
- Height: 35mm/1.38"

### Material Composition

- Latex & phthalate free
- 10mm foam
- Waterproof lining

### Order Information

| Product Code | Product Name                   | Qty    |
|--------------|--------------------------------|--------|
| CSM2608      | CS Bow & Wilson Comfort Covers | 6 Sets |

### Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

- If incorrectly applied to the Frame, this product can cause pressure marks or ulcers on the patient
- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624120244

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Ref: CSM2608 - TDS (Rev.03) 03/2026