



Technical Data Sheet

CSM2610 CS Jackson Standard Covers - 5 Piece



Description

- 2-Piece wedge cushion set
- For use on the Jackson Table
- Includes 1 chest cover, 2 hip covers and 2 thigh covers
- Helps to protect the chest and hip pads on the Jackson Table from wear and tear
- Aids infection control
- Light padding on the chest piece to reduce the risk of the patient experiencing marking or pressure ulcers

Application Area

- OR
- Anaesthesia

Measurements

- Length: 380mm/14.96"
- Width: 290mm/11.42"
- Height: 30mm/1.18"

Material Composition

- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2610	CS Jackson Standard Covers - 5 Piece	6 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- If incorrectly applied to the pads, the CS Jackson Table Compatible Covers can cause pressure marks or ulcers on the patient
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624120251

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Ref: CSM2610 - TDS (Rev.03) 03/2026