



Technical Data Sheet

CSM2714 CS Head Tray Attachment



Order Information

Product Code	Product Name	Qty
CSM2714	CS Head Tray Attachment	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Description

- Fits onto modular spinal tables to be used as a tray on which to place CS or other prone head supports
- Secure the tray to the side bars of the operating table and lock onto place with the quick-release locking clips
- Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)

Application Area

- OR
- Anaesthesia

Measurements

- Length: 500mm/19.69"
- Width: 300mm/11.81"
- Height: 60mm/2.36"

Material Composition

- ABS Plastic
- Latex & phthalate free

Cleaning

- Clean using non-alcohol hospital wipes

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Ensure the platform is securely attached to the table before use
- Check the platform for damage before and after each use. Do not use if damaged, contact supplier for advice
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624121098

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Australian Sponsor:

MDSS AU Pty Ltd, Ground Floor, Suite 55,
97-99 Bathurst Street, Sydney, NSW 2000, Australia

Ref: CSM2714 - TDS (Rev.03) 03/2026