



Technical Data Sheet

CSM2790 CS Patient Eye Protector



Description

- Protects the patients eyes against damages and dehydration during the surgery
- Transparent lenses allow an easy and comfortable control of the patient
- Self-adhesive and free of latex
- The flexible foam can be adapted for every face

Application Area

- OR
- Anaesthesia

Measurements

- Length: 210mm / 8.3"
- Width: 75mm / 3"
- Height: 10mm / 0.4"

Material Composition

- Material: Foam - made of PU
- Transparent shield - made of PE
- Glue: Medical pressure-sensitive adhesive (PSA)
- EO Sterelization Process, Non-allergenic

Order Information

Product Code	Product Name	Qty
CSM2790	CS Patient Eye Protector	25 units

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 25 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Avoid storage next to open doors and heaters

Contraindications

- Do not immerse in water
- Do not autoclave
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624124181

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark