



Technical Data Sheet

CSM4570 CS Pro-Axis Compatible Set



Description

- Single Use
- For use on the OSI ProAxis Table
- Helps to protect the chest and hip pads on the ProAxis Table from wear and tear
- Aids infection control
- Light padding on the chest piece to reduce the risk of the patient experiencing marking or pressure ulcers

Application Area

- OR
- Anaesthesia

Measurements

- Length: 640mm/25.2"
- Width: 640mm/25.2"
- Height: 2mm/0.08"

Material Composition

- Waterproof
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM4570	CS Pro-Axis Compatible Set	6 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

- If incorrectly applied to the pads, the CS ProAxis Table Compatible Covers can cause pressure marks or ulcers on the patient
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624120473

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark