



CARE SURGICAL


Product Code: CSM2820

CS Comfort Shin Supports (Pair)

⚠ Only to be used by trained clinicians. Ensure product is secure before positioning patient.

Instructions for use

- + Ensure feet do not touch operating table
- + Ankles should extend beyond Shin support



UK: Care Surgical Ltd
Unit 6 Ringtail Road,
Burscough,
L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



⚠ **SAFE DISPOSAL:** customers should adhere to any local laws and regulations to the safe disposal of medical devices.
⚠ **DO NOT USE** products that appear to be damaged or worn.
In an event of product damage or patient injury, please contact Care Surgical.



5 060624 120725

Ref: CSM2820 - IFU (Rev.09) 06/2025

