




CARE SURGICAL

Product Code: CSM2821

CS Knee Rings (Pair)

Instructions for use

- + Place the knees in the centre of the rings for minimum pressure



Only to be used by trained clinicians.

Ensure product is secure before positioning patient.


Cleaning instructions:



Clean using standard hospital wipes



Do not immerse in water



Do not autoclave



UK: Care Surgical Ltd
Unit 6 Ringtail Road,
Burscough,
L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



SAFE DISPOSAL: customers should adhere to any local laws and regulations to the safe disposal of medical devices.



DO NOT USE products that appear to be damaged or worn. In an event of product damage or patient injury, please contact Care Surgical.



Ref: CSM2821 - IFU (Rev.09) 03/2026