




Technical Data Sheet

CSM2411 CS Disposable Supine Head Support



Description

- Supports patient's heads in the supine position
- Single use - disposal after use with regular clinical waste
- Recommended for adults, but use will depend on patient's stature
- Reduces the incidence of pressure ulcers and marking

Application Area

- OR
- Anaesthesia

Measurements

- Length 220mm/9"
- Width 220mm/9"
- Height 70mm/2.8"

Material Composition

- Polyurethane foam core
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2411	CS Disposable Supine Head Support (Adult)	20 Units

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 20 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120343

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.03 28/01/2026