



Technical Data Sheet

CSM2520 CS ProneOne™ Head Support



Description

- Stand alone cushion to position the patients face in prone
- Universal face cut-out (suitable for patients from 8 years upwards)
- Optimal pressure distribution on forehead and cheeks
- Raises the patient's head to assist with achieving a natural, neutral neckline
- Mirror included to help monitor patient's eyes and ET tubes

Application Area

- OR
- Anaesthesia
- ICU

Measurements

- L: 282mm/11.10"
- W: 239mm/9.41"
- H: 138mm/5.43"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2520	CS ProneOne™ Head Support	8 Units

Packaging Information

- Each Unit is packaged in a recyclable polyethylene bag
- 8 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120039

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.03 05/03/2026