



Technical Data Sheet

CSM-2530 CS Prone Face Mask



Description

- Supports patients' faces in the prone position
- Compatible with the Allen C-Prone Head Positioner
- Distributes the pressure on the patient's cheeks, forehead and chin whilst in the prone position
- Reduces the incidence of pressure ulcers and marking

Application Area

- OR
- Anaesthesia

Measurements

- Length 330mm/12.99"
- Width 280mm/11.02"
- Height 50mm/1.97"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2530	CS Prone Face Mask	5 Units

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 5 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- POVL
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120572

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023