




# Technical Data Sheet

## CSM-2550V CS Prone Face Block 5"



### Description

- Soft foam block with slit/hole for ET tube on the right
- Firm base provides stability
- Soft foam top for patient comfort and protection
- Distributes the pressure on the patient's cheeks, forehead and chin while in the prone position
- Raises the patient's head to assist with achieving a natural, neutral neckline
- For patients aged 8 and over, dependent on stature
- For patient's weighing up to 200kg
- Remove from bag 1 (one) hour before use to allow the cushion to regain required shape and size
- Ensure product is fully inflated before use

### Application Area

- OR
- Anaesthesia

### Measurements

- L240 | W260 | H130 mm
- L9.4 | W10.2 | H5"

### Material Composition

- Polyurethane foam
- Latex & phthalate free

### Order Information

Product Code	Product Name	Qty
CSM-2550V	CS Prone Face Block 5"	14 units

### Packaging Information

- Each Unit is packaged in a recyclable polyethene bag (compressed)
- 14 units are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120091

### Labelling




Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 07/23