



# Technical Data Sheet

## CSM-2567 CS Adjustable Mirror Plus



### Description

- Mirror for use in the prone position
- The mirror enables the clinician to monitor the patient's eyes and tubing
- The cogs fitted to the mirror adjust the height or angle of the Prone Plus Helmet to achieve the safest natural neckline
- Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 220mm/8.66"
- Width 220mm/8.66"
- Height 4mm/0.16"

### Material Composition

- Made from ABS Plastic
- Latex & phthalate free
- Aluminium 6061

### Cleaning

- Clean the mirror with non-alcohol wipes or warm soapy water

### Order Information

Product Code	Product Name	Qty
CSM-2567	CS Adjustable Mirror Plus	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure ulcers
- POVL
- Neurological damage
- Venous pressure, which can lead to ischaemia
- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Not fully radiolucent- cog screws may be visible on x-rays
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120145

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023