




Technical Data Sheet

CSM2575 CS Arm Cushions (Pair)



Order Information

Product Code	Product Name	Qty
CSM2575	CS Arm Cushions (Pair)	20 Pairs

Packaging Information

- Each pair is packaged in a recyclable polyethene bag
- 20 pairs are packaged into a cardboard box for sale

Description

- Support patients' arms during procedure
- Reduces and distributes pressure on the patient's arms
- Reduces the risk of the patient experiencing marking or pressure ulcers
- Notches allow manipulation of the arm support to support and protect the patient's upper arm and forearm, or just the forearm

Application Area

- OR
- Anaesthesia

Measurements

- Length 600mm/23.62"
- Width 126mm/4.96"
- Height 65mm/2.56"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624-120183

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 08/2026