




Technical Data Sheet

CSM-2590 CS Adjustable Prone Leg Cushions



Description

- Supports patients' legs in the prone position
- 2 part product
- Part 1 supports the patient's knees (CSM 2590-1)
- Part 2 supports the patient's shins and allows their ankles/feet to find their own position (CSM 2590-2)
- Reduces the risk of pressure ulcers
- The 2 parts can be moved together or apart to suit the stature of the patient

Application Area

- OR
- Anaesthesia

Measurements

- Length 515mm/20.28"
- Width 500mm/19.69"
- Height 265mm/10.43"

Material Composition

- Polyurethane foam core
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2590	CS Adjustable Prone Leg Cushions	2 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 2 sets are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel
- Do not autoclave

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120213

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023