



Technical Data Sheet

CSM-2608 CS Bow & Wilson Comfort Covers



Description

- For use on the Allen Bow or Wilson Frame
- Disposable Drape with 10mm foam and waterproof lining
- Helps to protect the Allen Bow or Wilson Frame from wear and tear
- Aids infection control
- Light padding to reduce the risk of the patient experiencing marking or pressure ulcers

Application Area

- OR
- Anaesthesia

Measurements

- Length 760mm/29.92"
- Width 250mm/9.84"
- Height 35mm/1.38"

Material Composition

- Latex & phthalate free
- 10mm foam
- Waterproof lining

Order Information

Product Code	Product Name	Qty
CSM-2608	CS Bow & Wilson Comfort Covers	6 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- If incorrectly applied to the Frame, this product can cause pressure marks or ulcers on the patient
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120244

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023