



# Technical Data Sheet

## CSM-2812 CS Lateral Head Support (Adult)



### Description

- Supports patient's heads in the lateral position
- Reusable - clean using non-alcohol hospital wipes
- Recommended for adults, but use will depend on patient's stature
- Can be used for patients who are in the prone or lateral position
- Reduces the incidence of pressure ulcers and marking
- Patient can be repositioned without disconnecting their tubing, so PEEP pressure is maintained

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 208mm/8.19"
- Width 208mm/8.19"
- Height 51mm/2.01"

### Material Composition

- Polyurethane foam
- Latex & phthalate free

### Cleaning

- Clean using non-alcohol hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM-2812	CS Lateral Head Support (Adult)	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

### Contraindications

Incorrect positioning of the patient in the this product can increase the risk of:

- Pressure ulcers
- Neurological damage
- Venous pressure, which can lead to ischaemia
- Do not immerse in water
- Do not autoclave
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120404

### Labelling




Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023