



Technical Data Sheet

CSM2820 CS Shin Supports (Pair)



Description

- Reusable
- Large surface allows distribution of pressure on the patient's shins whilst in the prone position
- Reduces the incidence of pressure-related injury to the shins
- Non-slip base for secure positioning of the shins
- Helps to keep the patient's feet free of pressure

Application Area

- OR
- Anaesthesia

Measurements

- L240mm W240mm H180-250 mm
- L9.4" W9.4" H7.1-9.8"

Material Composition

- Polyurethane foam core
- Latex & phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2820	CS Shin Supports	1 Pair

Packaging Information

- Each pair is packaged in a recyclable polyethylene bag
- Each pair is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Nerve damage

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624-120725

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark