



Technical Data Sheet

CSM-2820 CS Shin Supports (Pair)



Description

- Reusable
- Large surface allows distribution of pressure on the patient's shins whilst in the prone position
- Reduces the incidence of pressure-related injury to the shins
- Non-slip base for secure positioning of the shins
- Helps to keep the patient's feet free of pressure

Application Area

- OR
- Anaesthesia

Measurements

- L240mm W240mm H180-250 mm
- L9.4" W9.4" H7.1-9.8"

Material Composition

- Polyurethane foam core
- Latex & phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM-2820	CS Shin Supports	1 Pair

Packaging Information

- Each pair is packaged in a recyclable polyethene bag
- Each pair is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Nerve damage

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120725

Labelling




Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 7/23