



Technical Data Sheet

CSM-2605 CS Wingset Comfort Covers



Order Information

Product Code	Product Name	Qty
CSM-2605	CS Wingset Comfort Covers	5 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 5 sets are packaged into a cardboard box for sale

Description

- For use with the Allen Wingset System
- Disposable Drapes with waterproof lining
- Helps to absorb excess moisture and protect the chest and hip pads on the Wingset System from wear and tear
- Aids infection control
- Light padding to reduce the risk of the patient experiencing marking or pressure ulcers

Application Area

- OR
- Anaesthesia

Measurements

- Length 234.5mm/9.23"
- Width 165mm/6.50"
- Height 10mm/0.39"

Material Composition

- White foam: Polyurethane
- Blue: Non-woven (waterproof): 55% Woodpulp 45% Polyester
- Latex & phthalate free

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- If incorrectly applied to the support pads, this product can cause pressure marks or ulcers on the patient
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120237

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023