




# Technical Data Sheet

## CSM2906 CS Tony Prone™



### Description

- Table top support for patients in the prone position for surgical and non-surgical procedures
- Fully radioluscent design- designed for imaging
- Reusable
- Distributes pressure along the torso the incidence of pressure related injury
- Recessed for abdominal decompression
- Comfortable for both conscious and unconscious patients

### Application Area

- OR
- Anaesthesia
- CT

### Measurements

- Length: 720mm/28.35"
- Width: 450mm/17.72"
- Height: 200mm/7.87"

### Material Composition

- Polyurethane foam core (coated)
- Latex & phthalate free

### Cleaning

- Clean with standard hospital wipes

### Order Information

| Product Code | Product Name   | Qty    |
|--------------|----------------|--------|
| CSM2906      | CS Tony Prone™ | 1 Unit |

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624122095

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark